

Patient Name  Date


Patient ID

Hospital Name  Doctor Name

Comments




### Gastroesophageal Junction Distensibility<sup>1</sup>



Diameter Area	
(mm)	(mm²)
5	20
6	28
7	39
8	50
9	64
10	79
11	95
12	113
13	133
14	154
15	177
16	201

#### EndoFLIP Instructions

- Step ① Just prior to inserting the catheter in the patient, place the balloon on a flat surface, and zero the pressure sensor with approx. 2-3ml in the balloon.
- Step ② Introduce the catheter to the **centre** of the GEJ.
- Step ③ Inflate the balloon to 30mL, and ensure an hourglass figure is obtained with the minimum diameter being to the middle part of the balloon. Note the balloon pressure (P) and minimum cross-sectional area (CSA min). Plot **CSA Min** versus **Pressure** on the graph above.
- Step ④ Repeat Step 3 for 20mL balloon volume.



<sup>1</sup>Increased Esophagogastric Junction Distensibility in GERD Patients Assessed with the Endoscopically Placed Functional Luminal Imaging Probe (Endoflip®)  
Monika A. Kwiatek, Ikuo Hirano, Peter J. Kahrlas, Barry P. McMahon, John E. Pandolfino Gastroenterology May 2009 (Vol. 136, Issue 5, Page A-16)

Patient Name  Date


Patient ID

Hospital Name  Doctor Name

Comments



### Gastroesophageal Junction Distensibility<sup>1</sup>



#### EndoFLIP Instructions

- Step ① Just prior to inserting the catheter in the patient, place the balloon on a flat surface, and zero the pressure sensor with approx. 2-3mL in the balloon.
- Step ② Introduce the catheter to the **centre** of the GEJ.
- Step ③ Inflate the balloon to 30mL, and ensure an hourglass figure is obtained with the minimum diameter being to the middle part of the balloon. Note the balloon pressure (P) and minimum cross-sectional area (CSA min). Plot **CSA Min** versus **Pressure** on the graph above.
- Step ④ Repeat Step 3 for 20mL balloon volume.



CSA min

<sup>1</sup>Increased Esophagogastric Junction Distensibility in GERD Patients Assessed with the Endoscopically Placed Functional Luminal Imaging Probe (Endoflip<sup>®</sup>)  
Monika A. Kwiatek, Ikuo Hirano, Peter J. Kahrlas, Barry P. McMahon, John E. Pandolfino Gastroenterology May 2009 (Vol. 136, Issue 5, Page A-16)