



NEW!

EndoFLIP®

Endolumenal Functional Lumen Imaging Probe



A revolutionary new technology used to measure the dimensions and function of the Lower Esophageal Sphincter



- **Testing of LES function**

EndoFLIP® is designed to test the function of the Lower Esophageal Sphincter by measuring its "distensibility".

- **LES distensibility**

The role of LES distensibility, a key determinant of opening dimensions, is mostly overlooked in current clinical assessment.

- **Real-time imaging of LES**

The EndoFLIP® system allows real-time imaging of LES intraluminal geometry during distention.

- **Dynamic quantitative measurement of LES**

EndoFLIP® allows dynamic quantitative measurement of LES distensibility while fitting in a practical clinical algorithm.

The EndoFLIP® Imaging System is currently CE-marked, but has not yet received FDA clearance. EndoFLIP® is not available for sales in the USA.

MMS
Medical Measurement Systems

P.O. Box 580
7500 AN Enschede
The Netherlands
T +31 - 53 - 480 37 00
F +31 - 53 - 480 37 01
E info@mmsinternational.com
I www.mmsinternational.com

M A K E S M E A S U R E M E N T S I M P L E